

# Repair of Distal Triceps Tendon Rupture with Human Acellular Dermal Matrix (ADM)

Andrea Atzei MD<sup>1</sup>, Giampietro Bertasi MD<sup>2</sup>

<sup>1</sup>University of Verona, Orthopedic Department, <sup>2</sup>San Pietro in Cariano, Verona

CASE STUDY

Distal triceps tendon ruptures can cause long-term disability for the patient, which can be further complicated by difficulties for the surgeon in developing a diagnosis as well as determining the severity of the injury.<sup>1,2</sup> Distal triceps tendon ruptures are often caused from a fall onto an outstretched hand or trauma to the posterior of the arm, with pulling or tearing at the osseous tendon insertion.<sup>1,3</sup> Other studies report intramuscular injury or injury at the myotendinous junction as another possible cause of the rupture.<sup>1,3</sup> Treatment of these ruptures typically includes surgery to reattach the ruptured triceps tendon to the olecranon of the elbow.<sup>4</sup> This is commonly achieved with the use of sutures that are passed through tunnels in the olecranon.<sup>4</sup> Surgeons often augment the tendon with an allograft in patients with chronic tears or undergoing revision surgery.<sup>5</sup>

An alternative treatment for ruptured distal triceps tendon is a matrix scaffold for new tissue generation, an acellular human dermal matrix (ADM) allograft as reviewed by Wainwright and Bury.<sup>6</sup> Decellularized human skin has been used for a variety of medical procedures, primarily involving wound healing, soft tissue reconstruction, and sports medicine applications.<sup>7-10</sup>

**The following case presentation involves treatment of a ruptured distal triceps tendon with this human decellularized dermis.**

## Patient

- 42 year-old male

## Diagnosis

- Failed repair of total full-thickness laceration of the distal tendon of the triceps (Fig. 1)
- Original failed repair had been attempted three months earlier

## Treatment

- The scar tissue was resected with a resulting 6 cm gap between the tendon stumps
- Two core sutures (Ethibond 2, Ethicon, Somerville, NJ, USA) were placed at both sides of the tendon (Fig. 2)
- One trimmed piece of 6cm x 7cm non-meshed decellularized dermis was weaved through the tendon mid-structure for augmentation and a second custom fit piece of 6cm x 7cm non-meshed decellularized dermis was laid over the repair site (Figs. 3,4)


- A small incision was made to let the exudate drain
- Elbow was splinted in 90° flexion for three weeks and progressive flexion was allowed to achieve full flexion in the following three weeks

## Conclusion


- Post-operative course was uneventful except for a prolonged swelling over olecranon bursa, healed spontaneously after three weeks
- By six months post-op, the patient regained 87% ROM and 70% strength
- Three month MRI demonstrated intact distal triceps brachialis tendon (Fig. 5)
- Excellent augmentation of the distal triceps brachialis tendon achieved

# Repair of Distal Triceps Tendon Rupture with Human Acellular Dermal Matrix (ADM)


CASE STUDY



**Figure 1.** Full-thickness laceration of the distal tendon of the triceps



**Figure 2.** Core sutures were placed at both sides of the tendon



**Figure 3.** Dermacell was weaved through the tendon mid-structure



**Figure 4.** A second custom-fit piece of Dermacell was laid over the repair site



**Figure 5.** Three month MRI demonstrates intact distal triceps brachialis tendons

## References

1. van Riet RP, Morrey BF, Ho E, O'Driscoll SW. Surgical treatment of distal triceps ruptures. *J Bone Joint Surg Am.* 2003;85A(10):1961-7
2. Rajasekhar C, Kakarlapudi TK, Bhamra MS. Avulsion of the triceps tendon. *Emerg Med J.* 2002;19:271-2
3. Sollender JL, Rayan GM, Barden GA. Triceps tendon rupture in weight lifters. *J Shoulder Elbow Surg.* 1998;7:151-3
4. Bava ED, Barber FA, Lunc ER. Clinical outcome after suture anchor repair for complete traumatic rupture of the distal triceps tendon. *Arthroscopy.* 2012;28(8):1058-63
5. Vidal AR, Allen A. Biceps tendon and triceps tendon injuries. *Clin Sports Med.* 2004;23(4):707-22
6. Wainwright, DJ, Bury, SB. 2011. Acellular Dermal Matrix in the Management of the Burn Patient. *Aesthet Surg J;* 31(7):135-235
7. Wong I, Burns J, Snyder S. 2010. Arthroscopic GraftJacket Repair of Rotator Cuff Tears. *J Shoulder Elbow Surg;* 19(2):104-109
8. Wilkins, R. 2010. Acellular Dermal Graft Augmentation in Quadriceps Tendon Rupture Repair. *Curr Orthop Pract,* 21(3):315-319
9. Lee, Daniel. 2007. Achilles Tendon Repair with Acellular Tissue Graft Augmentation in Neglected Ruptures. *J Foot Ankle Surg;* 46(6): 451-455
10. Sbitany, H., Sandeen, S., Amalfi, A., Davenport, M., Langstein, H. 2009. Acellular Dermis-Assisted Prosthetic Breast Reconstruction versus Complete Submuscular Coverage: A Head-to-Head Comparison of Outcomes. *Plast Reconstr Surg;* 124:1735-1740

LifeNet Health helps to save lives, restore health and give hope to thousands of patients each year. We are the world's most trusted provider of transplant solutions, from organ procurement to new innovations in bio-implant technologies and cellular therapies—a leader in the field of regenerative medicine, while always honoring the donors and healthcare professionals who allow the healing process.

## LifeNetHealth.org

The LifeNet Health logo is a registered trademark of LifeNet Health.  
©2019 LifeNet Health, Virginia Beach, VA. All rights reserved.  
68-20-097.03