

# Combined Aorto-Caval Reconstruction Using Two Cryopreserved Aortoiliac Allografts Following Retroperitoneal Tumor Resection

Case performed by: Amit B. Karmur, D.O.<sup>1</sup>; Jeremy Blumberg, M.D.<sup>2</sup>; Nina M. Bowens, M.D.<sup>1</sup>

CASE STUDY

<sup>1</sup>Division of Vascular and Endovascular Surgery, Harbor-UCLA Medical Center, Torrance, CA

<sup>2</sup>Division of Urology, Harbor-UCLA Medical Center, Torrance, CA

## Abstract

Oncologic resection of retroperitoneal tumors often requires complex vascular reconstruction. Autologous and/or synthetic grafts are most commonly used; however, they can be associated with significant morbidity. This case describes the successful outcomes and vascular reconstruction of the infrarenal aorta and inferior vena cava using two cryopreserved bifurcated aortoiliac (AI) allografts following resection of a large retroperitoneal tumor with vascular involvement.

## Introduction

Large retroperitoneal tumors can be very challenging to manage, with surgical resection as the primary treatment option and often the only hope for patient survival. However, the involvement of major organs and blood vessels has been considered a contraindication for oncological resection due to high surgical risk and the potential for poor clinical outcomes.<sup>1-4</sup> Nonetheless, reconstruction of major blood vessels, including the aorta and inferior vena cava (IVC), following *en bloc* resection of large tumor masses can be achieved with positive clinical outcomes in carefully selected patients.<sup>2-8</sup>

Autologous veins and synthetic grafts are two of the most commonly used graft options for blood vessel reconstruction. However, these options can be associated with significant morbidity, particularly in cases when there is a risk of contamination and infection.<sup>1, 6, 9-13</sup> Cryopreserved human allografts are an alternative option for *in situ* arterial reconstruction, which offer the advantage of avoiding a second surgical site while being resistant to infection.<sup>14, 15</sup> Aortoiliac (AI) allografts, in particular, offer a more complex anatomy through the presence of branch vessels, making them an advantageous choice for reconstructing major vessels like the aorta and IVC. This case describes the successful combined reconstruction of both the infrarenal aorta and IVC using two cryopreserved AI allografts in the

setting of retroperitoneal tumor resection. The patient consented to publication of this case report and images.

## Case Report

A 20-year-old male patient presented with persistent and worsening abdominal pain, unintentional weight loss, vomiting, and early satiety over a period of three months, after having had no previous health concerns. A testicular mass and abdominal distension were noted upon physical examination. Computed tomography (CT) scans revealed a 15 cm necrotic, multi-lobulated retroperitoneal mass with near-circumferential involvement of the aorta and iliac bifurcation and compression of the IVC (Figure 1).



Figure 1. Initial presentation computed tomography (CT) image showing a 15 cm multi-lobulated tumor mass with near-circumferential involvement of the infrarenal aorta and IVC. The IVC was compressed, and there was thrombus noted in both the IVC and right common iliac vein.

Additionally, there was thrombus noted in both the IVC and right common iliac vein. The patient was diagnosed with metastatic non-seminomatous testicular cancer. Following a full course of neoadjuvant chemotherapy, the tumor regressed to 6 cm, but still showed persistent involvement of the aortic bifurcation and the IVC (Figure 2).

# Combined Aorto-Caval Reconstruction Using Two Cryopreserved Aortoiliac Allografts Following Retroperitoneal Tumor Resection

## CASE STUDY



Figure 2. Preoperative CT image following neoadjuvant chemotherapy. The tumor mass reduced to 6 cm, but still showed involvement of the aorta and IVC.

With the persistent tumor mass involvement of both the aorta and IVC, reconstruction of both blood vessels was planned in conjunction with a retroperitoneal lymph node dissection. Cryopreserved bifurcated AI allografts (LifeNet Health®, Virginia Beach, VA) were chosen for the reconstruction of both the aorta and the IVC given the possibility of greater tumor involvement than perceived on CT, potentially including the venous bifurcation as well as the bowels. In preparation for the vessel reconstructions, measurements of the infrarenal aorta and IVC were made in the axial and coronal planes on the preoperative CT scan (Figure 2). The infrarenal aorta diameter was 13 mm, and the right and left common iliac arteries were 9 mm each. Because the tumor's significant effect on the IVC made planning for its reconstruction more challenging, diameter measurements were taken from the cephalad portion of the mass and each common iliac vein, which were 19 mm and 11 mm, respectively. Based on these measurements, the diameters of the AI allografts used for the aorta and IVC reconstructions were 15 mm and 21 mm, respectively, when distended under pressure (125 mmHg). Distended diameter is used to emulate the diameter of the recipient's vessel under systemic pressure *in vivo*, enabling more accurate caliber matching. Using a distended AI graft diameter of 21 mm for reconstruction of the IVC, which experiences a much lower systemic pressure, would allow for a higher velocity within the venous lumen *in vivo* and decrease the risk of graft thrombosis.

The patient underwent laparotomy and was systemically heparinized prior to vessel clamping and resection. *En bloc* resection was then performed including the residual tumor with nodes along with the infrarenal aorta with bilateral common iliac arteries and the IVC with bilateral iliac veins. Reconstruction of the IVC and bilateral iliac veins was performed first (Figure 3A), followed by the infrarenal aorta and bilateral iliac arteries (Figure 3B).




Figure 3. Intraoperative images taken after the reconstruction of the A) IVC and B) infrarenal aorta, each using a cryopreserved AI allograft. The distended diameter of the AI allograft used to reconstruct the IVC was 21 mm, and that for reconstructing the aorta was 15 mm to match the patient's native vessel diameters.

Each reconstruction was performed in an end-to-end configuration for both proximal and distal anastomoses. Immediately postoperative, the patient had palpable peripheral pulses as well as Doppler signals in the proximal vena cava. The patient did not experience any complications and was discharged from the hospital within one week. One-month postoperative imaging revealed patent reconstructions, as well as complete tumor mass resection (Figure 4).

# Combined Aorto-Caval Reconstruction Using Two Cryopreserved Aortoiliac Allografts Following Retroperitoneal Tumor Resection

## CASE STUDY



Figure 4. Postoperative CT image taken after 1 month demonstrating complete tumor resection including the involved blood vessels, as well as patent reconstructions of both the IVC and aorta.

## Discussion

Major vascular involvement is traditionally considered a contraindication for oncological resection. However, with precise preoperative planning in carefully selected patients, resection and vascular reconstruction can be achieved. This case demonstrates successful combined aorto-caval reconstruction using two cryopreserved AI allografts in the setting of retroperitoneal tumor and lymph node resection.

Various graft options exist for the reconstruction of major blood vessels. Autologous vein grafts and synthetic grafts are two of the most commonly used; however, each is associated with limitations. Autologous veins, often considered the preferred option<sup>11</sup>, require harvesting from a second surgical site, which can be a source of pain and possible infection, and there may be limited availability in some patients due to poor vascularity. Additionally, the potential for caliber mismatch between the autologous vein and the native vessel is a concern.<sup>6, 9, 10</sup> Synthetic grafts are another common choice; however, the use of synthetic grafts should be avoided when contamination is a risk.<sup>1, 13</sup> Cryopreserved allografts present an alternative option, offering the advantages of avoiding a second surgical site and being resistant to infection.<sup>14, 15</sup>

As with autologous grafts, ensuring appropriate caliber matching between the allograft and native vessel is essential for achieving a successful and patent reconstruction. Singular allograft vessels often require surgical modifications in order to match the patient's vessel diameter and anatomy, which can increase operative time. Cryopreserved AI allografts, through the presence of the branch vessels, can form direct anastomoses with nearby vessels, providing conduits for reconstruction of complex vascular anatomy.<sup>15</sup> Moreover, because AI allografts distend when exposed to arterial pressure (the extent of which can be variable between grafts) measuring the pressurized distended versus flaccid diameter of the allograft reduces the possibility of caliber mismatch upon implantation. The AI allografts used in this case were measured under distended pressure (125 mmHg) by LifeNet Health, allowing for accurate caliber-matching during the preoperative planning and eliminated the need for surgical modifications. This resulted in the AI implants being well-matched anatomically with the patient's native vessels, leading to patent reconstructions and successful clinical outcomes, overall.

## Conclusion

This case study describes the positive clinical outcomes in a challenging case of a large retroperitoneal tumor with involvement of major blood vessels, and the successful combined reconstruction of the infrarenal aorta and IVC using two cryopreserved AI allografts. While not generalizable, the results from this case support the use of cryopreserved AI allografts for major vessel reconstruction following *en bloc* oncological resection.

# Combined Aorto-Caval Reconstruction Using Two Cryopreserved Aortoiliac Allografts Following Retroperitoneal Tumor Resection

## CASE STUDY

## References

1. Ghosh J, Bhowmick A, and Baguneid M. Oncovascular surgery. *Eur J Surg Oncol* 2011;37(12):1017-24.
2. Newton DH, Chen SL, Fraker DL, Guzzo TJ, Roses RE, Foley PJ, et al. PC214. Techniques and Outcomes of Concomitant Aortic and Caval Resection and Reconstruction for Cancer. *Journal of Vascular Surgery* 2019;69(6):e263-e264.
3. Carpenter SG, Stone WM, Bower TC, Fowl RJ, and Money SR. Surgical management of tumors invading the aorta and major arterial structures. *Ann Vasc Surg* 2011;25(8):1026-35.
4. Bower TC, Nagorney DM, Cherry KJ, Toomey BJ, Hallett JW, Panneton JM, et al. Replacement of the inferior vena cava for malignancy: An update. *J Vasc Surg* 2000;31:270-81.
5. Bertrand MM, Carrere S, Delmond L, Mehta S, Rouanet P, Canaud L, et al. Oncovascular compartmental resection for retroperitoneal soft tissue sarcoma with vascular involvement. *J Vasc Surg* 2016;64(4):1033-41.
6. Fernandez-Alonso L, Fernandez-Alonso S, and Aguilar EM. Iliac Artery Reconstruction With the Superficial Femoral Vein After Major Oncologic Resection. *Vasc Endovascular Surg* 2017;51(5):350-351.
7. Schneider JR, Sener SF, and Barrera E, Jr. Combined replacement of infrarenal aorta and inferior vena cava after en bloc resection of retroperitoneal extraosseous osteosarcoma. *J Vasc Surg* 2008;48(2):478-9.
8. Schwarzbach MH, Hormann Y, Hinz U, Leowardi C, Bockler D, Mechtersheimer G, et al. Clinical results of surgery for retroperitoneal sarcoma with major blood vessel involvement. *J Vasc Surg* 2006;44(1):46-55.
9. Siddiqi S. Saphenous vein harvest wound complications: risk factors, identification, prevention, and management. *Chronic Wound Care Management and Research* 2016;Volume 3:147-156.
10. Tennyson C, Young CP, and Scarci M. Is it safe to perform endoscopic vein harvest? *Interact Cardiovasc Thorac Surg* 2010;10(4):625-9.
11. O'Connor S, Andrew P, Batt M, and Becquemin JP. A systematic review and meta-analysis of treatments for aortic graft infection. *J Vasc Surg* 2006;44(1):38-45.
12. Clagett GP, Valentine RJ, and Hagino RT. Autogenous aortoiliac/femoral reconstruction from superficial femoral-popliteal veins: Feasibility and durability. *J Vasc Surg* 1997;25:255-70.
13. Wilson WR, Bower TC, Creager MA, Amin-Hanjani S, O'Gara PT, Lockhart PB, et al. Vascular Graft Infections, Mycotic Aneurysms, and Endovascular Infections: A Scientific Statement From the American Heart Association. *Circulation* 2016;134(20):e412-e460.
14. Brown KE, Heyer K, Rodriguez H, Eskandari MK, Pearce WH, and Morasch MD. Arterial reconstruction with cryopreserved human allografts in the setting of infection: A single-center experience with midterm follow-up. *J Vasc Surg* 2009;49(3):660-6.
15. Harlander-Locke MP, Harmon LK, Lawrence PF, Oderich GS, McCready RA, Morasch MD, et al. The use of cryopreserved aortoiliac allograft for aortic reconstruction in the United States. *J Vasc Surg* 2014;59(3):669-74.

Results from case studies are not predictive of results in other cases. Results in other cases may vary.

LifeNet Health helps to save lives, restore health and give hope to thousands of patients each year. We are the world's most trusted provider of transplant solutions, from organ procurement to new innovations in bio-implant technologies and cellular therapies—a leader in the field of regenerative medicine, while always honoring the donors and healthcare professionals that allow the healing process.

## LifeNetHealth.org

LifeNet Health and the LifeNet Health logo are registered trademarks of LifeNet Health.

©2021 LifeNet Health, Virginia Beach, VA. All rights reserved.  
68-20-324.00